

INTRODUCING...



THE **23** MINUTE HOUR

SIMPLY AMAZING.

With a 10 inch chamber and closed door drying, the Zipclave sterilizes and dries instruments in approximately 23* minutes or less.



ZIPCLAVE

Summit Dental Systems introduces another reliable workhorse in the field of dentistry: the Zipclave Sterilizer. Once set, the Zipclave will automatically go from heating cycle, to sterilization cycle, to the closed-door vent/dry cycle in approximately 18-23* minutes. With no cool down period between cycles, the Zipclave allows instant back-to-back cycles, becoming more efficient after every use.



Using steam sterilization, the Zipclave works faster due to the instant-dry system, a process that takes less than 2 minutes to produce dry material at the end of each cycle. With a 3 year limited parts warranty, the Zipclave provides the peace of mind one needs to run their practice efficiently.

*So what are you going to do
with your extra time?*

* in a warm to warm cycle

CLOSED-DOOR DRYING USING INSTANT DRY

Does not expose the instruments to airborne pathogens during the drying process, which takes less than 2 minutes

QUICK STERILIZATION AND DRYING CYCLES

Complete sterilization and drying cycles in approximately 18-23 minutes

- No need for an additional/secondary sterilizer
- Higher patient turnover rate per instrument used
- Less instruments need to be purchased/used

NO MICROELECTRONICS

- No cool down period between cycles
- Speeds up subsequent sterilizations

SIMPLE DESIGN FEATURES

Uncomplicated gauges and cycle setting parameters, with pre-programmed sterilization cycles for effortless setting (determined by load type)

LITTLE TO NO CLEANING

Water is purged after every cycle, which does not allow biofilm to build up in the water reservoir

SAFETY SYSTEM

A double locking door mechanism guarantees the door cannot be opened while the chamber is pressurized

STERILIZATION CYCLE DURATION

6 (unwrapped), 10 or 20 (wrapped) minutes - (as selected cycle)

DIMENSIONS

24" x 17" x 15" (608 x 431 x 380 mm)

CHAMBER VOLUME

23 liters

CHAMBER DIAMETER

10 in (254mm)

CHAMBER DEPTH

18 in (457.2mm)

WATER CONSUMPTION PER CYCLE

8fl oz (237ml)

WATER RESERVOIR CAPACITY

1.1 gallons (4 liters) for the first supply

PRESSURE STERILIZATION

6 or 10 minutes: 30PSI (2.10 kgf/cm²) / 20 minutes: 20PSI (1.40 kgf/cm²)

OPERATING TEMPERATURE

6 or 10 minutes: 275° F (135° C) / 20 minutes: 257° F (125° C)

NET/GROSS WEIGHT

86 lb (39 Kg) / 94 lb (42.5 Kg)

VOLTAGE/FREQUENCY

120 VAC/60 Hz

POWER

1600 VA

ZIPCLAVE

So what are you going to do with your extra time?



800.275.3368 | Tel: 954.730.3636 | Fax: 954.730.3602
6277 NW 28th Way | Fort Lauderdale, FL 33309

www.summitdental.com